


- Augments and restores midfacial deficiency
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Easy to trim with scissors or scalpel
- Screws or sutures can be placed in any desired location
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4156	Small - Left	28mm	26mm	4.5mm
4157	Small - Right	28mm	26mm	4.5mm
4158	Large - Left	30mm	28mm	7mm
4159	Large - Right	30mm	28mm	7mm

