



Surgical Implants

Designed with Michael A. Burnstine, M.D.

- “Excellent addition to the management of patients with symptomatic paralytic ectropion” - Michael A. Burnstine, M.D.
- Thin, tissue-like profile for greater eyelid stability
- Fully porous to promote tissue ingrowth and facilitate implant stabilization
- Placed in Lamellar Space inferior to the tarsus under no tension
- Available with strip to provide additional support
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



Scan for Video Demonstration



CAT#	Description	A	B	C
4436	Regular	29.5mm	9mm	0.45mm
4437	Regular + Superior Strip	29.5mm	9mm	0.45mm