


Surgical Implants

- Restores volume of the orbit following enucleation and evisceration procedures
- Fully porous to promote tissue growth
- Wrapping of implant is optional, as extraocular muscles may be sutured directly to sphere
- Introducer provided sterile with each implant
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



Introducer



Supreme suturability without the need of wrapping.

CAT#	Diameter
4028	14mm
4029	16mm
4030	18mm
4031	19mm
4032	20mm
4033	21mm
4034	22mm
4035	23mm

