



### Surgical Implants

- Su-Por Biomaterial Sheets provide the surgeon with excellent options for craniofacial reconstruction and augmentation
- The individually packaged, sterile implants provide “off-the-shelf” availability, and may save the time and expense of harvesting graft material
- Su-Por Sheets are available in a variety of sizes and in thicknesses ranging from 0.25mm to 3.0mm
- Feathering the edge of the sheet allows for a more smooth transition from the implant to the adjoining bone
- Retained Shape Memory™ provides malleability by hand
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	A	B	C
4001	38mm	50mm	0.25mm
4002	50mm	76mm	0.25mm
4003	38mm	50mm	0.3mm
4004	50mm	76mm	0.35mm
4005	30mm	50mm	0.40mm
4006	38mm	50mm	0.45mm
4007	50mm	76mm	0.45mm
4008	38mm	50mm	0.85mm
4009	50mm	76mm	0.85mm
4012	38mm	50mm	1.5mm
4013	50mm	76mm	1.5mm
4015	38mm	50mm	3.0mm

