



## Surgical Implants

- Provides anterior projection
- Designed to be trimmed to meet the needs of the patient
- Fully porous material for tissue ingrowth
- Implant bends and holds desired shape without heating
- Fixation screws or sutures can be placed in any desired location
- Easy to trim with scissors, fluted burr, or scalpel
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	Description	A	B	C
4064	Left	43mm	18mm	3.2mm
4065	Right	43mm	18mm	3.2mm