


Surgical Implants

- Provides an alternative to tissue grafting
- No donor site morbidity
- May carve implants in the sterile field to obtain contours individualized for the surgical situation
- Allowing the implant to soak in a hot, sterile saline bath relaxes the memory of the implant, enabling easier modification into desired shape
- Individually packaged, MR Safe, Sterile, 5 Year Shelf Life



CAT#	A	B	C
4019	13mm	38mm	3mm
4020	25mm	50mm	3mm
4021	38mm	63mm	3mm
4022	13mm	38mm	6mm
4023	25mm	50mm	6mm
4024	38mm	63mm	6mm
4025	13mm	38mm	9.5mm
4026	25mm	50mm	9.5mm
4027	38mm	63mm	9.5mm

